# bewell<sup>®</sup> connect

✓RoHS

CE


# MyScale

BW-SC1

By Visionmed

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BW-SC1

By Visionmed

This product complies with the Reach standard (RoHS).

## TABLE OF CONTENTS

1. NOTES ON SAFETY .....	23
2. INFORMATION .....	23
2-1. ABOUT BODY FAT .....	24
2-2. ABOUT BODY WATER .....	24
2-3. ABOUT MUSCLE MASS .....	24
2-4. ABOUT BONE MASS .....	25
2-5. ABOUT BMI - BODY MASS INDEX .....	25
2-6. ABOUT DAILY CALORIE REQUIREMENTS .....	26
3. USE .....	26
3-1. DESCRIPTION .....	26
3-2. PRECAUTIONS BEFORE USE .....	26
3-3. SCALES USER GUIDE .....	27
3-4. DOWNLOAD APPLICATION .....	31
3-5. DATA TRANSFER .....	31
3-6. CHANGING THE BATTERIES .....	31
4. HOW TO MAINTAIN AND STORE THE UNIT .....	31
5. HOW TO PREVENT A MALFUNCTION .....	32
6. SPECIFICATIONS .....	32
7. TROUBLESHOOTING .....	34

Log onto our website to find out how to install and use the BewellConnect® application

[www.bewell-connect.com/install](http://www.bewell-connect.com/install)

*The manufacturer keeps the right to modify without any preliminary opinion technical specifications of the product.*

Dear customer,

Thank you for buying the connected and intelligent BewellConnect® body analyser scales. We hope you enjoy using them and we recommend you read the instructions carefully.

## 1. NOTES ON SAFETY

- Use this appliance only for the use for which it is designed as described in this manual.
- This appliance may be used for personal use at home. Do not use the unit for other purpose.
- Use this appliance in an ambient temperature range of 0 to 40°C.
- Do not expose this appliance to extreme temperature conditions > 60°C ou < -20°C.
- Do not use this appliance at a relative humidity of more than 85%.
- This appliance must always be placed in a clean and dry place.
- Do not expose this appliance to sunlight or water.
- Do not expose this appliance to electric shocks.
- Do not use this appliance outdoors.
- Never drop the appliance.
- Follow the maintenance instructions specified in this manual.
- Do not attempt to open the device. In case of problems, contact your dealer.
- Do not leave this appliance within the reach of children.
- The measurement results are given for information. If in doubt about your results, please contact your doctor.
- Discontinue use of the device in case of anomalies or malfunction.
- Do not use any battery other than the mentioned batteries. Do not attempt to recharge non-rechargeable batteries. Do not throw them into a fire.
- Remove the batteries if the appliance is not used for a long period.
- This appliance is not designed for use by people (including children) with reduced physical, emotional or mental abilities or those lacking experience or knowledge unless used in the presence of someone responsible for their security, under surveillance or following instructions about how to use the appliance.

### Important

*Observe the polarity of the batteries. Incorrect polarity can cause damage and compromise the warranty of your appliance.*

## 2. INFORMATION

The human body is almost entirely made up of the following: bones, muscles, organs, fat and water.

## 2-1. ABOUT BODY FAT

Body fat protects organs, controls body temperature and contains vitamins and energy. It is not unhealthy but excess body fat can cause health problems. The body fat percentage provides an insight into a person's physical condition.

The perfect body fat percentage (also known as adiposity) depends on age and gender.

STANDARD BODY FAT PERCENTAGE (for reference only)								
Age	Female				Male			
	Slim	Normal	Over-weight	Obese	Slim	Normal	Over-weight	Obese
15 – 39	14 – 21	22 – 27	28 – 32	33 – 40	10 – 18	19 – 24	25 – 30	31 – 40
40 – 55	15 – 22	23 – 28	29 – 33	34 – 40	11 – 19	20 – 25	26 – 31	32 – 40
56 – 70	16 – 23	24 – 29	30 – 34	35 – 40	12 – 20	21 – 26	27 – 32	33 – 40

*Please speak to your doctor if you have any questions.*

## 2-2. ABOUT BODY WATER

Water is an essential part of the body. The amount of water the body contains reflects how healthy a person is. Water plays a vital role in almost all the organs in the human body.

It gives cells form, helps maintain body temperature, transports nutrients throughout the body, hydrates skin and improves digestion. It helps remove waste from the body, protect vital organs and lubricate joints.

The perfect body water percentage depends on age and gender.

STANDARD BODY WATER PERCENTAGE (for reference only)		
Age	Female	Male
10-15	57-67%	58-72%
16-30	47-57%	53-67%
31-60	42-52%	47-61%
61-80	37-47%	42-56%

*We recommend you speak to a doctor if you have any questions. The results may be unreliable if the subject is dehydrated.*

## 2-3. ABOUT MUSCLE MASS

Muscle plays a vital role as it acts as your body's engine by consuming energy (calories). Muscle mass is based on the weight of the body's muscle (skeletal muscle, smooth muscle and water contained in muscle).

### STANDARD MUSCLE MASS PERCENTAGE (for reference only)

Female	Male
34%	40%

## 2-4. ABOUT BONE MASS

Bone mass is the weight of the skeleton's bones. It accounts for approximately 15% of total weight. The percentage gradually decreases from the age of 50. We recommend you speak to a doctor in the event of rapid bone loss.

### STANDARD BONE MASS IN KG (for reference only)

Weight	Female
< 50	1,95
50 - 75	2,40
> 75	2,95
Weight	Male
< 65	2,66
65 - 95	3,29
> 95	3,69

## 2-5. ABOUT BMI - BODY MASS INDEX

BMI measures body fat and is an indicator of potential health risks based on weight/height ratio.

### BMI INDICATOR (for reference only)

BMI	Diagnostic
< 16.5	Severely underweight
16.5 - 18.5	Underweight
18.5 - 25	Normal
25 - 30	Overweight
30 - 35	Moderately obese
35 - 40	Severely obese
> 40	Morbidly obese

## 2-6. ABOUT DAILY CALORIE REQUIREMENTS

All food provides calories and all calories provide energy. The human body uses the energy to live, breathe and carry out every day activities. Men consume more energy than women. In general, men require 5-10% more energy than women with the exception of pregnant or breastfeeding women whose daily calorie intake ensures mother and child are in good health. Based on your age, height and gender, the BewellConnect® scales will tell you how many calories you need every day to achieve a healthy weight.

## 3. USE

### 3-1. DESCRIPTION

SET : On/Off  
in Weight and Calculator mode  
Profile configuration  
Confirmation button

Adjustment of input  
data value

### 3-2. PRECAUTIONS BEFORE USE

- The appliance is not suitable for pregnant women, people with a pacemaker, people on dialysis, people with an implant (e.g. prosthetics, contraceptive implants, pin, breathing apparatus) as the appliance may malfunction. Please speak to your doctor if you have any questions.
- The appliance is not suitable for children aged under 10 to calculate body fat, water mass and bone mass. Do not leave children unattended when using the appliance.
- The appliance is not suitable for people aged over 100 to calculate body fat, water mass and bone mass.
- The appliance is not suitable for professional athletes or bodybuilders (intensive bodybuilding). Results displayed may be inaccurate due to physiological variations, muscle density and athletes' hydration levels. Please speak to your doctor if you have any questions.
- Use the scales on a sturdy, flat non-slip surface (e.g. not on a rug).
- Only stand on the scales when 0:0 appears on the screen.
- You must be barefoot (no shoes or socks) or else the result may not be displayed or may be unreliable.
- Weight and body mass varies throughout the day. We recommend you weigh yourself before meals for comparable results. It is advisable to analyse weight and body mass at the same time every day, preferably on an empty stomach.
- The results displayed by this appliance are for reference only and not provided for medical purposes. If your results are higher or lower than normal, please ask your doctor for advice.

## 3-3. SCALES USER GUIDE

The BewellConnect® scales have two modes:

- Weight Mode
- Impedancemeter (Weight and Calculation Mode : body fat, muscle, water and bone mass and daily calorie requirements).

- |                        |                          |                        |
|------------------------|--------------------------|------------------------|
| 1 Age                  | 2 Muscle mass icon       | 3 Bone mass icon       |
| 4 Gender (man/woman)   | 5 Profile number         | 6 Weight status        |
| 7 Height unit          | 8 Weight unit            | 9 Hydration percentage |
| 10 Body fat percentage | 11 Body mass index (BMI) | 12 Calories            |

**Overload warning :** Remove items from scales immediately if *Err* appears on the screen. This may cause permanent damage to the appliance.

### Weight unit selection

Turn appliance on. When 0:0 appears, press the *UNIT* button on the back of the appliance to selection Kilos, Pounds or Stones.

### Weight mode

Stand on the scales. The scales turn on. 0:0 appears on the screen followed by your weight. The numbers flash then stabilise and stay on.

If you want to weigh yourself without selecting a profile, the scales will only display weight. If you create your profile, you will be able to view the health indicators measured by the scales displayed twice in this order: weight, FAT% (body fat percentage), TBW% (total water mass percentage), Muscle mass, Bone mass, BMI (body mass index), Kcal (daily calories requirements to maintain weight).

## Personal data set-up

Press the **SET**  button to set up your user profile. The screen displays the flashing profile number.


Press the  /  buttons to select a user profile (1 - 10). Press **SET**  to confirm your selection.

Press the  /  buttons to select your gender (male/female). Press **SET**  to confirm your selection.

Press the  /  buttons to enter your age. Press **SET**  to confirm your selection.

Press the  /  buttons to enter your height (in centimetres if measurement unit is kilos (kg) / feet and inches if measurement unit is pounds (lb) and stones (st)). Press **SET**  to confirm your selection.

## Impedancemeter mode (Weight and calculation)


To calculate your body fat, turn on the scales by pressing **SET**  or by placing your foot on it. The last profile used is displayed on the screen with the age and height of the user.

Choose the user profile using the  /  p buttons. The age and height of the selected profile appears on the screen followed by 0:0.

Stand on the balance to weigh yourself.

For accurate results, you must place your feet on the scales' sensors and remain stationary whilst being weighed.

Once complete, the results appear on the screen and are repeated.

WEIGHT


FAT% : BODY FAT PERCENTAGE

TBW% : WATER MASS PERCENTAGE

### 3-4. DOWNLOAD APPLICATION

Download the BewellConnect® app on your mobile or tablet:

- App store or Google Play

- scan the QR Code on the side of the box  
Then click on the MyScale icon.

### 3-5. DATA TRANSFER

Connect your appliance (mobile, tablet) using Bluetooth (Bluetooth 4.0 smart Low Energy connectivity suitable for iPhone 4S) to track your weight and measurements: data transfers automatically with the BewellConnect® app. Simply add the measurement to a profile. View your profile history on the app.

### 3-6. CHANGING THE BATTERIES

The battery compartment is at the back of the appliance.

- Insert 4 new AAA batteries in the scales before use.

- When the icon appears on the LCD screen, the batteries are low or uncharged. Replace old batteries and make sure the positive and negative ends are facing the correct direction as an error may cause damage or affect your appliance's warranty.

Never insert rechargeable batteries. Use single-use batteries. Do not mix old and new batteries.

Disposal of empty battery to the authorized collecting party subject to the regulation of each individual territory. If the unit will not be used for a long period, please take out all batteries to avoid liquid leakage damaging the unit.

### 4. HOW TO MAINTAIN AND STORE THE UNIT

- Keep scales dry. Water exposure may make them slippery and cause an accident.

- Store your scales horizontally.

- Store your scales in a dry place away from humidity.

- Use a soft and dry cloth to clean appliance. If necessary, slightly dampen cloth with tap water.

- Do not use alcohol, benzene, diluent or other chemical products to clean appliance.

- This product is not waterproof. Avoid contact with liquid.

- Please leave at room temperature for 30 minutes before use if appliance has been exposed to temperatures under 0°C otherwise results may be unreliable.

- Please remove batteries (they may damage appliance) if appliance is unused for a long time. Make sure the positive and negative ends are facing the correct direction. Do not mix old and new batteries. Do not use rechargeable batteries. Do not dispose of used batteries in normal household waste.

## 5. HOW TO PREVENT A MALFUNCTION

- Do not drop the device and avoid sudden jars or shocks.
- Do not insert other objects into any holes.
- Do not attempt to disassemble the unit.
- Do not put the unit in water.
- The monitor might not meet its performance specifications or cause safety hazard if stored or used outside the specified temperature and humidity ranges in specifications.

## 6. SPECIFICATIONS

Product name	BewellConnect® MyScale		
Model	BW-SC1		
Data	Units	Minimum	Maximum
Weight measuring range	kg	6	150
BMI range		2.5	90
Body fat range	%	5	50
Water mass range	%	15	85
Muscle mass range	%	5	60
Bone mass range	kg	0	10
Calories range	Kcal	385	3999
Graduation accuracy			
Accurate to 50 kg	kg	0.3	
Accurate to 100 Kg	kg	0.5	
Accurate to 150 Kg	kg	0.7	
Percentage accuracy			
Body fat	%	5	
Water mass	%	5	
Muscle mass	%	5	
Bone mass accuracy	kg	0.2	
Calories accuracy	Kcal	115	
Weight units	kg / lb / st		

Profiles	10 user profiles Age : 5-99 years old Height : 100-200 cm
Dimensions	330 x 330 x 22mm
Weight	2kg
Working Temp. & humidity	Temp : 0 - 40°C - Humd : 85%
Storage Temp. & humidity	Temp : -10 - 60°C - Humd : 85 %RH
Bluetooth technology	4.0 smart Low Energy
Bluetooth distance	Max. 10m
Consumption	≤20mA
Power source	4 x AAA batteries (provided)
Battery life	3-4 months (based on 4 measurements per day) Low battery indicator
Automatic shut down	±7 seconds
Supplied with	A connected and intelligent body analyser, user manual, quick installation guide, 4 AAA batteries

## 7. TROUBLESHOOTING

If you have one of the following problems, please refer to this breakdown service guide to help resolve the problem. If the problem persists, please contact our customer service.

PROBLEM	POSSIBLE CAUSE	SOLUTION
Nothing appears when you press SET  or when you press the scales	No battery inserted	<i>Insert batteries.</i>
The  icon appears	Used batteries	<i>Insert new batteries.</i>
	Overload warning	<i>The scales' maximum capacity is 150kg. Remove weight from scales immediately as this may permanently damage the appliance.</i>
	The body fat or muscle mass percentage cannot be found.	<i>Communication error during measurement. Please ensure you are barefoot whilst being weighed.</i>
	Less than 5% body fat	
	More than 50% body fat	

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**EN:** Visiomed\* will repair or replace this product free of charge in the case of defective parts or manufacturing defects, in accordance with the conditions mentioned below as follows:

**DURATION:** 24 MONTHS RETURN TO WORKSHOP

**LIMITS AND EXCLUSIONS:** This guarantee concerns only the original final purchaser. A purchase invoice, or another proof of purchase, with this guarantee card will be required to obtain an after-sales service, in accordance with this guarantee. This guarantee card will not be extended to another person only the original final purchaser. This guarantee becomes void if the serial numbers on the product are modified, replaced, illegible, absent, or if repair has been carried out by a service not approved, including the user.

This guarantee covers only the defects of the material or parts, occurring during normal use of the product. It does not cover the damage caused during the transport of the apparatus, causes due to repairs being carried out by the distributor, by any modifications undertaken, any connection of equipment not approved by Visiomed\*, or causes contrary to those written in the user manual or notice. Moreover, the present guarantee does not cover damage due to falls, bad handling, bad installations, damage by fire, floods, lightning, or any other natural disaster. This guarantee does not cover the packing of the material, the accessories, the defects caused by commercial exposure of the product, show room, sale space, demonstration etc... Normal maintenance, cleaning and the replacement of parts where wear is normal, are not covered by the terms of this guarantee. Visiomed\* and its representatives and agents will not in any case be held responsible for any damage and consecutive damages due to the mishandling of this product. This guarantee is the only valid one at Visiomed\*, any other guarantee (commercial guarantee) except this one will not be taken into account.

**IMPORTANT:** During the guarantee period if you are dissatisfied with the repairs of this product, please contact the Visiomed\* customer service.

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**EN:** Disposal of Waste Equipment by Users in Private Households in the European Union.

The symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

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© 06/2015 BW-SC1
Visiomed SAS France